




gmmedical
Product Catalogue



GM Medical's extensive range of equipment continues the tradition of industry leadership.


With the recent inclusion into the ParagonCare Group GM Medical has expanded the range of equipment to offer significantly more products into the healthcare environment.

GM Medical has developed a comprehensive equipment range that combines advanced functionality and contemporary design with quality and reliability.

Using a variety of finishes and colour combinations the product range offers unlimited applications throughout the healthcare environment, from the stringent infection control requirements of the operating theatre and CSSD right through to the Aged Care facilities and Consulting Rooms with colour coordinated environments.

CONTENTS

Instrument Trolleys	01
Equipment Trolleys	04
Spectra Trolleys	08
Sixty Forty Materials Management	09
CSSD Equipment	10
Operating Theatre	12
General Ward	16
Maternity	24
Linen and Supply	26
Food Services	31
Consulting Rooms	32
Lighting	37
Product Guide	38



GM Medical instrument trolley design methods are class leaders in providing high quality carts for the surgical and clinical environments.

With a focus on durability features such as double folded safety edges, under shelf bracing and castor mount inserts as standard. These ensure safety and stability at all times.



Shelf Rails

12mm stainless steel rails
Fully welded stainless steel shelf
Double safety folds to all edges

Under Shelf Bracing

Full under bracing providing maximum stability and strength

Drawer Runners

Full extension, easy glide metal runners with drawer retention system

Stainless Steel Construction

Satin finish - 304 grade stainless steel throughout

Shelf - Stainless Steel

Double safety folds to all edges

Castors

75mm ball bearing fork castor
Soft grey rubber wheel
Expanding adaptor for easy removal

INSTRUMENT TROLLEYS

NO RAIL



AX 001
Instrument Trolley
490 x 490 x 900mm



AX 003
Instrument Trolley
600 x 490 x 900mm



AX 006
Instrument Trolley
750 x 490 x 900mm

RAIL 3 SIDES



AX 033
Instrument Trolley
490 x 490 x 900mm



AX 036
Instrument Trolley
600 x 490 x 900mm



AX 039
Instrument Trolley
750 x 490 x 900mm

ONE DRAWER



AX 054
Instrument Trolley
490 x 490 x 900mm



AX 057
Instrument Trolley
600 x 490 x 900mm



AX 060
Instrument Trolley
750 x 490 x 900mm

TWO DRAWER



AX 066
Instrument Trolley
490 x 490 x 900mm



AX 069
Instrument Trolley
600 x 490 x 900mm



AX 072
Instrument Trolley
750 x 490 x 900mm



AX 009
Instrument Trolley
 900 x 490 x 900mm
 Other models available
 AX 0015
 1000 x 600 x 900mm
 AX 0018
 1200 x 600 x 900mm



AX 042
Instrument Trolley
 900 x 490 x 900mm
 Other models available
 AX 0048
 1000 x 600 x 900mm
 AX 0051
 1200 x 600 x 900mm



AX 063
Instrument Trolley
 900 x 490 x 900mm



AX 075
Instrument Trolley
 900 x 490 x 900mm



AX 076
Instrument Trolley
 490 x 490 x 900mm
 Stainless steel construction
 3 Drawers - 125mm deep
 125mm castors

Other models available
 AX 077
 600 x 490 x 900mm
 3 Drawers - 125mm deep
 AX 0780
 750 x 490 x 900mm
 3 Drawers - 125mm deep
 AX 079
 900 x 490 x 900mm
 3 Drawers - 125mm deep



AX 080
Instrument Trolley
 490 x 490 x 900mm
 Stainless steel construction
 4 Drawers - 125mm deep
 125mm castors

Other models available
 AX 081
 600 x 490 x 900mm
 4 Drawers - 125mm deep
 AX 082
 750 x 490 x 900mm
 4 Drawers - 125mm deep
 AX 083
 900 x 490 x 900mm
 4 Drawers - 125mm deep




AX 073
Instrument Trolley
 900 x 490 x 900mm
 Stainless steel construction
 2 half drawers - 125mm deep
 75mm castors



AX 107
Resuscitation Trolley
 750 x 490 x 900mm
 Stainless construction
 4 half drawers - 125mm deep
 125mm castors

EQUIPMENT TROLLEYS



The range of Equipment Trolleys is extensive and all trolleys are designed featuring a solid tubular frame which ensures rigidity at all times. This allows the drawers to be mounted securely and to operate efficiently. All drawer runner systems are full ball bearing extension type featuring a push lock mechanism to prevent the drawer from accidentally opening during transport. As with the instrument trolleys the

shelves feature double fold and under shelf bracing for a stable work platform. Our range of Equipment trolleys are available in all stainless steel or can be colour coordinated with powder coated drawers and panels.

All trolleys are fitted with 125mm precision bearing heavy duty castors with optional corner buffers available.



AX 108
Resuscitation Trolley - Standard
 750 x 490 x 900mm
 Stainless construction
 4 Half drawers - 125mm deep
 1 Full width drawer - 125mm deep
 125mm castors



AX 110
Resuscitation Trolley
 750 x 490 x 900mm
 Stainless construction
 6 Half drawers - 125mm deep
 1 Full width drawer - 250mm deep
 125mm castors



AX 078
Multi Drawer Trolley
 900 x 490 x 900
 Stainless steel construction
 4 Drawers - 100mm deep
 2 Drawers - 150mm deep
 125mm castors



AX 104
Anaesthetic Trolley
 600 x 490 x 1000mm
 Stainless construction
 2 Drawers - 125mm deep
 2 Drawers - 250mm deep
 125mm castors



AX 102
Anaesthetic Trolley
 600 x 490 x 900mm
 Stainless construction
 3 Drawers - 125mm deep
 1 Drawer - 250mm deep
 125mm castors



AX 105
Anaesthetic Trolley
 900 x 490 x 900mm
 Stainless construction
 6 Full width drawers
 125mm castors

AX 114 (not pictured)
Anaesthetic Trolley
 600 x 490 x 900mm
 Stainless construction
 5 Full width drawers -
 125mm deep
 125mm castors

EQUIPMENT TROLLEYS



AX 117
Emergency Equipment Trolley
750 x 490 x 900mm
Stainless steel frame
Powder coated panels and drawers
4 Drawers - 105mm deep
1 Drawer - 225mm deep
125mm castors



AX 115
Emergency Equipment Trolley
750 x 490 x 900mm
Stainless steel frame
Powder coated panels and drawers
5 Drawers - 125mm deep
125mm castors



AX 145
IV Set-up Trolley
900 x 490 x 1122mm
5 Drawers - 125mm deep
2 Half drawers - 250mm deep
125mm castors



AX 119
Emergency Equipment Trolley
750 x 490 x 1120mm
Stainless steel frame
Powder coated panels and drawers
3 Drawers - 125mm deep
2 Drawers - 250mm deep
125mm castors

MEDICATION TROLLEYS



AX 676
Medication Trolley - 3 Drawer
 Stainless construction
 600 x 490 x 1100mm
 Lift up shelf
 3 Drawers - 275mm deep
 125mm castors



AX 677
Medication Trolley - 1 Drawer
 Powder coated finish
 600 x 490 x 900mm
 1 Webster drawer
 275mm deep - lockable
 Lift up shelf
 75mm castors

TROLLEY ACCESSORIES



AX 120
Drawer insert containers - 12 units



AX 738
Self closing drawer runner



AX 132
Storage panel - 4 sml / 3 med / 1 lrg



AX 126
Drop down work surface



AX 123
Vertical oxy bottle holder and clamp



AX 129
IV pole and clamp



AX 135
4 x pin power board



AX 143
Push handle
 to suit 400mm wide trolley



AX 142
Waste bin and bracket

SPECTRA TROLLEYS

The Spectra range is a trolley system that combines mobility and drawer flexibility into a lightweight colour coded package that has applications throughout all healthcare environments.

The unique drawer systems is based on a removable 60 x 40 tray system that offers 3 different depths and the option of integrated drawer dividers.

This allows flexibility with the trolley applications, particularly in the ability to restock the carts with transfer tray systems.

The shelves are one piece moulded polymer which are easy to clean and eliminates infection control issues due to the rounded internal mould design, whilst integrating a raised lip to three sides to prevent items falling off.



AX S0001
2 Shelf Trolley
 700 x 500 x 1000mm
 Work surface area 620 x 460mm
 125mm castors

 Push handle optional on all



AX S0010
1 Drawer Trolley
 700 x 500 x 1000mm
 AX M0601 - 1 x 100mm deep tray
 Interlock dividers optional
 Work surface area 620 x 460mm
 125mm castors



AX S0020
2 Drawer Trolley
 700 x 500 x 1000mm
 AX M0601 - 2 x 100mm deep trays
 Interlock dividers optional
 Work surface area 620 x 460mm
 125mm castors



AX S0030
Multi Drawer Trolley
 700 x 500 x 1000mm
 AX M0605 - 1 x 50mm deep tray
 AX M0601 - 2 x 100mm deep trays
 AX M0602 - 2 x 200mm deep trays
 Interlock dividers optional
 Work surface area 620 x 460mm
 125mm castors

COLOUR OPTIONS

RED

YELLOW

GREEN

BLUE



60 x 40 Tray System



50mm Trays

AX M0605 50mm deep tray
 AX M1665 Interlock divider 600mm
 AX M1645 Interlock divider 400mm



100mm Trays

AX M0601 100mm deep tray
 AX M1661 Interlock divider 600mm
 AX M1641 Interlock divider 400mm



200mm Trays

AX M0602 200mm deep tray
 AX M1662 Interlock divider 600mm
 AX M1642 Interlock divider 400mm



The 60 x 40 materials management system is recognised throughout the world as the industry leader in hospital and healthcare storage. GM Medical brings these advanced storage solutions to Australia.

The 60 x 40 concept is based on a tray size of 600 x 400mm and a choice of tray depths of 50, 100, 200mm. This system allows versatility and flexible

configuration, as the trays are interchangeable with the range of GM cabinets, shelving systems and carts.

Choose your Configuration All GM Medical 60 x 40 trolleys provide you the opportunity to configure them to your requirements. Mix and match using our suggestions pictured as your starting point.



AX M2000
Mobile Storage Unit
 900 x 650 1000mm
 2 x 50 / 6 x 100 / 2 x 200mm deep trays
 125mm castors - 2 brake
 Push handle



AX M1400
Mobile Storage Unit
 665 x 460 x 975mm
 4 x 100 / 1 x 200mm deep trays
 100mm castors - 2 brake
 Push handle



AX M2500
Mobile Storage Unit
 Stainless Steel Construction
 600 x 460 x 900mm
 2 x 120mm deep drawers
 1 x 200mm deep tray
 125mm castors



AX M8000
Mobile Storage Unit
 680 x 450 x 1950mm
 1 x 50 / 4 x 100mm deep trays
 2 x 200mm deep trays
 125mm castors - 2 brake



AX M8500
Static Storage Unit
 680 x 450 x 1950mm
 1 x 50 / 4 x 100mm deep trays
 2 x 200mm deep trays
 Adjustable feet



Tray Stop System
 Prevents accidental removal of trays.
 Multi position holes enable flexible configuration.

CSSD EQUIPMENT

Our range of CSSD workstations and associated equipment are designed with flexibility in mind.

A large range of accessories, options and sizes are available and allows the workstation to be configured to meet a multitude of applications, optimising workflow, room layout and departmental efficiency.

The simple push button electric height adjustment ensures that the optimum working height is achieved to ensure correct operator's posture, providing comfort, and minimising occupational health and safety issues.

The workstations are available in six standard sizes, however should there be a requirement for

a particular size, we are able to custom manufacture.

Construction of all stainless steel, featuring a reinforced heavy duty top surface, ensures infection control issues are minimised. Stable base design ensures stability and provides clear access for operators legs and feet.

PICTURED WITH ACCESSORIES



Accessory Shelf

Constructed of chrome plated wire. 175mm wide. May be positioned at different heights if required.

Option - Shelf depth 300mm

Additional shelves may be added

AX 0672 1800 x 175mm

AX 0677 1800 x 300mm

AX 0673 1500 x 175mm

AX 0678 1500 x 300mm

AX 0674 1200 x 175mm

AX 0679 1200 x 300mm

AX 0685

Sterile Wrap Rack

This rack system can be fitted either side of the workstation and is removable and telescopic.

Option - The wrap holders have the option of being telescopic to accommodate wider rolls of sterile wrap.

AX 0687

Keyboard Shelf

Pull out keyboard shelf on quality roller drawer runners. May be positioned anywhere along the length of the workstation.

AX 0686

CPU Holder

Robust construction. Can be positioned anywhere beneath the workstation surface.

Designed to for standard CPU units.

Option - custom size CPU holder

Additional Accessories (not pictured)

AX 0689 Lower Storage Shelf

AX 0690 Power Board Outlet 6 plug

AX 0691 Steri-Peel Holder

Size Options

AX 0666 1800 x 760mm

AX 0669 1800 x 900mm

AX 0667 1500 x 760mm

AX 0670 1500 x 900mm

AX 0678 1200 x 760mm

AX 0671 1200 x 900mm

BASIC UNIT



Electric Height Controls

Simple push button controls conveniently located at the corner of the workstation.

Dual column electric height adjustment.

Standard height range 800-1100mm.

The height range is 300mm.

Option - Battery back up system



AX 0688
Preparation Table Adjustable Height
 1800 x 900mm / height range 830 - 1130mm
 Electrically operated height adjustment
 Stainless steel work surface
 Storage shelves / battery back up
 100mm braking castors



AX 692
Case Cart
 900 x 600 x 1040mm
 Stainless steel construction
 Fold out hinged doors to 270°
 Internal adjustable wire shelves
 Top rail push handle incorporated
 Full bumper surround
 200mm castors - 2 brake



AX 278
List Cart
 1200 x 600 x 1430mm
 Stainless steel construction
 2 hinged doors / internal adjustable shelves
 200mm castors / full bumper surround
 In addition to the standard sizes offered,
 custom sizes are available for both List
 and Case Carts.

Cart Features



Door System
 The doors swing open to a full 270° and ensures full unobstructed access. Flush fitting, heavy duty door latches are fitted. The doors have bracing rods welded to the door perimeter to ensure rigidity and durability. The door hinges are all stainless steel, powder coated white.

Shelf System
 Each unit is supplied with an adjustable shelf system that allows the stainless steel wire shelf to be positioned as required. Additional shelves are available.



AX 0684
Transportation Trolley
 1200 x 600 x 900mm
 Stainless steel construction
 3 shelves with rails to three sides
 125mm castor




AX 0681
Kinguard Drying Rack
 1200 x 500 x 1500mm
 Stainless steel construction
 6 rails / 125mm castors



AX 0680
Packing Table
 Stainless steel construction
 1 drawer - 125mm deep
 Pull out keyboard shelf

OPERATING THEATRE



The GM Medical range of operating theatre equipment has a strong emphasis on being manufactured from stainless steel minimising infection control issues. All products in this range are fully welded and highly polished, therefore easy to clean and maintain.

This diverse range of products covers all clinical and surgical equipment requirements.



AX 106
Airway Trolley
 750 x 490 x 900mm
 Stainless construction
 1 Drawer - 125mm deep
 Clear Acrylic top



AX 4990
Plaster Trolley
 900 x 490 x 900mm
 Stainless construction
 2 Drawers - 125mm deep
 Drop down shelf / bucket hook



AX 252
Mayo Table
 Height 820 - 1120mm
 Tray size 560 x 410mm
 Stainless construction
 2 Leg flat base
 75mm castors



AX 253
Mayo Table
 Height 820 - 1120mm
 Tray size 560 x 410mm
 Stainless construction
 2 Leg slide base
 75mm castors



AX 254
Mayo Table
 Height 820 - 1120mm
 Tray size 560 x 410mm
 Stainless construction
 3 Leg base
 75mm castors



AX 256
Mayo Table
 Height 820 - 1120mm
 Tray size 560 x 410mm
 Stainless construction
 4 Leg base
 50mm castors



AX 021
Nesting Tables - Set of Three
 AX 024
Nesting Table - Large
 1000 x 600 x 900mm
 AX 027
Nesting Table - Medium
 940 x 510 x 880mm
 AX 030
Nesting Table - Small
 840 x 490 x 840mm



AX 111
Arm Table - Fixed Height
 750 x 490 x 900mm
 Stainless steel construction
 Safety folds to all edges
 Reinforced shelf
 2 x 75mm castors



AX 113
Arm Table - Adjustable Height
 750 x 490mm
 Height range 740 - 1140
 Stainless steel construction
 Reinforced shelf
 2 x 75mm castors
 Radiolucent top



AX 260
Lithotomy Table
 750 x 490mm Adjustable height
 4 - 75mm castors

OPERATING THEATRE



AX 273
Surgeons Stool



AX 277
Theatre Stool

AX 273

Surgeons Stool - Hydraulic

- Special operators stool for surgery
- Height adjustment by maintenance free hydraulic pump lift
- Safety base made of stainless steel with 5 antistatic twin castors
- Chromium plate foot ring
- Stroke 190mm
- Seat height range 530 - 720
- Back rest adjustment range 70mm
- Upper part swivelling & lockable
- Comfortable antistatic black upholstery
- Pair of upholstered OP - arm supports
- Base locking device

AX 277

Theatre Stool

- Height adjustment by safety gas spring controlled by stainless steel wobble plate
- Multi purpose stool for surgery, anaesthesia and outpatient department
- Safety base made of stainless steel with 5 twin castors for hard floors
- Round seat with antistatic black upholstery
- Seat height range 500 - 680mm
- Base locking device all round operatable by step disc above base



AX 262
Surgeon Stool Hand Operated
Height 530 - 710mm
No backrest



AX 265
Surgeon Stool Hand Operated
Height 530 - 710mm
Backrest



AX 267
Surgeon Stool Foot Operated
Height 530 - 710mm
No backrest



AX 270
Surgeon Stool Foot Operated
Height range 530 - 710
Backrest



AX 2710
Saddle Stool
Hand operated height adjustment
Height range 530 - 730mm



AX 391
One Step Stool with hand rail
Aluminium tread plate
700 x 300 x 200mm
Handle height; 1000mm



AX 387
One Step Stool Aluminium Tread
400 x 300 x 200mm
Stainless steel frame
Aluminium tread plate

AX 405
Two Step Stool Aluminium Tread
Step Height: 200 / 400mm
Stainless steel frame
Aluminium tread plate

AX 390
One Step Stool Aluminium Tread
700 x 300 x 200



AX 214
Double Bowl Stand
 Height 820mm
 300mm (dia) Stainless steel bowl



AX 211
Single Bowl Stand
 Height 820mm
 300mm (dia) Stainless steel bowl



AX 248
Single Bowl Stand
 Height 515mm
 300mm (dia) Stainless steel bowl



AX 240
Kick Bucket Stand
 Bucket and stand combination
 Stainless steel construction
 Separate components
 AX 241 Kick Bucket Stand
 AX 244 Kick Bucket - 11.4 litre
 AX 247 Kick Bucket Lid



AX 701



AX 702

AX 701
ICU Work Station
 Writing surface: 750 x 620mm
 Height range: 700 - 1000mm
 Powder coated
 Tilting writing surface

AX 702
ICU Work Station
 1500 x 600 x 1000mm
 Powder coated
 Tilting writing surface
 Pull out keyboard shelf
 Storage basket
 125mm castors - 2 brake
 Option - adjustable height
 writing surface



AX 841
Modular AV Equipment Trolley
 (monitor not included)
 680 x 550 x 1500mm
 Height to top shelf - 1100mm
 Pull out keyboard shelf
 100mm castors




AX 840
Modular Equipment Trolley
 680 x 550 x 1000mm
 Height to top shelf - 1100mm
 1 Drawer
 100mm castors



AX 636
Theatre Supply Storage Cart
 6 Swing out baskets - 2 fixed
 Laminated top
 4 - 100mm twin castors

GENERAL WARD



This range of equipment is engineered for strength and durability to withstand continuous use year after year. A blend of different materials and finishes are used throughout the range.

All products are fully welded for maximum rigidity and are supported by a full 2 year warranty. A choice of all stainless steel or colour coordinated powder coating is available.



AX 215
Pogo IV Pole / Stand - Heavy Base
 Adjust height by twist mechanism
 Height range 1350 - 2100mm
 5 Leg base - 75mm castors



Lift lock adjustment



AX 2380
Provita Smart IV Stand
 One hand operation
 Height range 1300 - 2100
 Stable five star weighted base
 4 hook IV holder
 Colour coordinated
 Aqua seal / 80mm castors



One hand adjustment




AX 231
Storage Basket and Clamp



AX 232
Vertical Oxy Holder and Clamp




AX 233
IV Handle and Clamp



AX 217
IV Pole / Stand
 Height range 1350 - 2050mm
 Stainless steel construction
 5 Leg weighted base
 50mm twin castors



AX 219
IV Pole / Stand - Heavy Base
 Height range 1350 - 2050mm
 Stainless steel construction
 5 Leg weighted base
 50mm twin castors



AX 216
IV Pole / Stand Polymer Base
 Height range 1350 - 2100mm
 5 Leg base



AX 220
IV Transfusion Pump Stand
 Height range 1350 - 2100mm
 Offset pole for stability
 Push handle

GENERAL WARD



AX 528
Oxygen Bottle Trolley G Size
 Epoxy coated - 150mm castors
 AX 531

Oxygen Bottle Trolley G Size
 Stainless steel - 150mm castors



AX 546
Oxygen Bottle Trolley D Size
 Epoxy coated - 150mm castors
 AX 549

Oxygen Bottle Trolley D Size
 Stainless steel - 150mm castors



AX 537
Oxygen Bottle Trolley C Size
 Epoxy coated - 150mm castors
 AX 540

Oxygen Bottle Trolley C Size
 Stainless steel - 150mm castors



AX 524
Multi Stack Gas Cylinder Trolley
 Capacity 9 x C size cylinders
 AX 525

Multi Stack Gas Cylinder Trolley
 Capacity 12 x C size cylinders



AX 524SS (pictured above)
Multi Stack Gas Cylinder Trolley
 Capacity 6 x D size / 6 x C size cylinders
 Stainless steel construction
 AX 525SS

Multi Stack Gas Cylinder Trolley
 Capacity 12 x C size cylinders
 Stainless steel construction



AX 519
Garbage Trolley-Foot Operated
 Foot operated lid action
 Stainless construction - PVC lid



AX 177
Lead Apron Rack
 700 x 600 x 1420mm
 Stainless steel construction
 12 apron hooks / corner buffers
 125mm castors



AX 175
Garment Rack
 1600 x 600 x 1600mm
 Stainless steel construction
 AX 178 Coat Hangers



AX 990
Charts Trolley - 12 Charts
 465 x 530 x 900mm
 Stainless steel construction
 125mm castors



AX 993
Charts Trolley - 24 Charts
 805 x 530 x 900mm
 Stainless steel construction
 125mm castors



AX 655
Files Trolley
 750 x 490 x 900mm
 Stainless steel construction
 1 Drawer - 125mm deep
 Removable clear dividers - Lift up side shelf
 125mm castors



AX 996
Chart Holder A4
 To suit Charts Trolleys
 310 x 265mm

AX 994
Charts Trolley - 36 Charts
 1170 x 530 x 900mm



AX 554
Pan and Bottle Rack - 4
 Capacity 4 pans / 5 bottles
 Stainless construction
 Includes driptray and back plate



AX 552
Pan and Bottle Rack - 6
 Capacity 6 pans / 7 bottles
 Stainless construction
 Includes driptray and back plate



AX 501
Bed Cradle
 Tubular frame
 Stainless steel construction



AX 534
Gas Bottle Holder
 130 x 130 x 350
 Hook over / to suit C size
 Sheet Metal



AX 513
Urinal Bottle Holder
 Stainless Steel
 AX 510
Urinal Bottle Holder
 Nylon Coated



AX 522
Garbage Wall Mount Bracket
 310 x 190
 Stainless steel
 Pre drilled mounting holes

AX 543
Gas Bottle Holder
 150 dia x 330
 Hook over / to suit C size
 12mm dia rod to hook into
 trolley or bed

GENERAL WARD

AX 432

Shower Chair

The shower chair offers a three position, drop down armrest, a non-slip retracting footrest fitted with an anti-tilt bracket and comes with a seamless, soft padded seat and backrest. The shower chair's slim line design is made of stainless steel and is fitted with water resistant precision bearing castors which offers excellent manoeuvrability.

Features

- Seat width 440mm
- Seamless padded seat and backrest
- Bed pan holder
- Padded arm rest
- Locking castors
- Stainless steel construction

Options

- Oxygen cylinder holder and clamp
- IV pole and clamp
- Hinged seat
- Safety arms - wrap around

Optional Accessories

(selection shown right)

AX 436 Hinged Padded Seat

AX 438 Padded Enclosed Seat

AX 441 Padded Open Front Seat

AX 444 Padded Backrest

AX 462 Oxygen Cylinder Holder & Clamp

AX 460 IV Pole and Clamp

AX 455 Safety Arms - Wrap Around

AX 447 Leg Extension



AX 432SA

Shower Chair

Swing back arms





AX 434
Shower Chair
 Seat width 440mm / drop arms
 Large rear 150mm castors



AX 433
Shower Chair
 Seat width 600mm / drop arms
 Weight capacity 175kg



AX 435
Shower Chair - Bariatric
 Seat width 660mm / swing back arms
 Weight capacity 350kg / heavy duty frame
 Reinforced seat / backrest

ACCESSORIES



AX 460
Shower Chair
 IV Pole and Clamp



AX 462
Shower Chair
 Vertical Oxy Holder and Clamp



AX 455
Shower Chair
 Safety Arms - Wrap Around



AX 447
Shower Chair
 Leg Extension

GENERAL WARD



AX 708

Fixed Height Patient Trolley

2000 x 685 x 750mm

Stainless steel construction

Removable side rails and end rails

Mediflex mattress

Storage basket / includes oxygen bottle holder

200mm castors / 3 brake / 1 steer



AX 704

Hi Low Mortuary / Concealment Trolley

2200 x 650mm / height range 450 - 1450mm

Stainless steel construction

Hi low function electrically operated by hand controller / battery back up for transportation

Removable side rails and end rails

200mm castors / 3 brake / 1 steer



AX 707

Mortuary / Concealment Trolley Cover

Document Holders

GM document holders offer a variety of configurations for the effective storage of files, binders and documents.

The range offers the following options:

- Landscape positioning
- Portrait positioning
- Multiple document holders
- Magnetic mounting option
- Clear or dark tinted options
- Unbreakable polycarbonate construction for some models



AX 9144

Document Holder

Clear - 260 x 254 x 102mm

Tilt Bins

Modular units with stackable design

- Manufactured using High Impact Polystyrene Plastic (HIPS)
- Transparent bins keep supplies visible
- Supplied with product bin labels for easy identification



AX 9142

Document Holder

Clear polycarbonate - 254 x 267 x 51mm



AX 9140

Document Holder

Clear polycarbonate - 445 x 165 x 76mm



AX 9136

Document Holder - Set of 3

Clear polycarbonate - 368 x 165 x 76mm



AX 9146

Document Holder

Clear - 375 x 178 x 102mm



AX 9122

Tilt Bins - 6 Units



AX 9124

Tilt Bins - 5 Units



AX 9126

Tilt Bins - 4 Units




AX 9128

Tilt Bins - 3 Units



MATERNITY



The GM range of Maternity products provide for varying conditions of infants, including the ability to bathe. Our range features a twin bassinet which assists in the bonding of the newborns. With Occupational Health and Safety issue being at the forefront we offer 2 types of variable height bassinets units

that offer simple electrically controlled height adjustments with a convenient hand control. This ensures safe, secure handling and treatment of the infant.



AX 609

Single Bassinet Trolley - Adjustable Height

790 x 420mm / height 800 - 1100mm
Electrically operated by hand controller / battery back up
Includes acrylic tub / mattress and 1 x swing out basket



AX 610

Babies Bath - Adjustable Height

790 x 420mm / height 800 - 1100mm
Stainless steel construction
Electrically operated by hand controller / battery back up
Includes acrylic tub with waste outlet connection and hose
1 x swing out basket
AX 611 - Acrylic tub (only) with waste outlet connection and hose



AX 582

Bassinet Trolley - Single

800 x 415 x 910mm
Includes acrylic tub and mattress
2 x white storage trays
Powder coated frame
AX 607 - Stainless steel frame



AX 597

Bassinet Trolley - Twin

800 x 610 x 910mm
Includes acrylic tub and mattress
2 x white storage trays
Powder coated frame
AX 608 - Stainless steel frame




AX 612

Bassinet Trolley - Cabinet

800 x 415 x 910mm
Cabinet with internal adjustable shelf
Includes acrylic tub and mattress

LINEN AND SUPPLY



This extensive range of equipment offers handling and transport solutions for linen and other supplies within the healthcare environment. Engineered to withstand heavy loads and frequent use, we have designed these products to be ergonomic

and easy to use. All equipment is fully welded and has a sturdy frame design for maximum durability and is constructed from stainless steel or durable powder coated steel.



AX 300 DOUBLE LINEN SKIP
WITH OPTIONAL HALF HEIGHT SHELF



Foot operated lid



AX 319

Soft Close Lid System - Optional

This unique mechanism provides a soft closing action which ensures smooth and quiet operation. This feature can be fitted to all linen skips with PVC lids.



AX 315

Swivel Hooks - Optional

4 hooks per linen bag ring available on all linen skips in the range



AX 291

Single Linen Skip - foot operated lid

380 x 540 x 850mm
Stainless steel construction
Foot operated PVC lid
75mm castors



AX 300

Double Linen Skip - foot operated lid

780 x 540 x 850mm
Stainless steel construction
Foot operated PVC lid
75mm castors



AX 307

Triple Linen Skip - foot operated lid

1170 x 540 x 850mm
Stainless steel construction
Foot operated PVC lid
75mm castors

Colour Identification Strips

Apply to PVC lid for colour coding (all models)

- AX 318R Red
- AX 318B Blue
- AX 318G Green
- AX 318Y Yellow

LINEN AND SUPPLY



AX 303
Triple Linen Skip - no lid
 1170 x 540 x 850mm
 Stainless steel construction
 75mm castors



AX 315
Swivel Hooks - Optional
 4 hooks per linen bag ring
 available on all linen skips
 in the range



AX 294
Double Linen Skip - no lid
 780 x 540 x 850mm
 Stainless steel construction
 75mm castors



AX 286
Single Linen Skip - no lid
 380 x 540 x 850mm
 Stainless steel construction
 75mm castors



AX 320
**Single Linen Skip
 Half height shelf
 Optional**

POWDER COATED



AX 307 PC
Triple Linen Skip - foot operated
 1170 x 540 x 850mm
 Powder coated finish
 Foot operated PVC lid
 75mm castors



AX 300 PC
Double Linen Skip - foot operated
 780 x 540 x 850mm
 Powder coated finish
 Foot operated PVC lid
 75mm castors



AX 291 PC
Single Linen Skip - foot operated
 380 x 540 x 850mm
 Powder coated finish
 Foot operated PVC lid
 75mm castors



AX 733
Moist Laundry Trolley
 Including spring loaded base
 450 litre capacity



AX 729
Moist Laundry Trolley
 Including spring loaded base
 320 litre capacity



AX 091
Utility Trolley
 720 x 450 x 910mm
 Stainless steel frame
 2 lift out tubs - 200mm deep
 125mm castors



AX 726
Wet and Dry Laundry Trolley
 1000 x 650 x 820mm
 Including spring loaded base



AX 731
Spring loaded base only



AX 094
Storage Trolley
 700 x 450 x 1300mm
 Stainless steel frame
 4 slide out tubs - 200mm deep
 125mm castors



AX 634
Bulk Collection Trolley
 Removable tray
 Stainless steel construction
 125mm castors



AX 0730
Soiled Linen Trolley
 1150 x 575 x 1240mm
 Powder coated finish
 Gas assisted lid (as shown)
 Full bumper surround / 200mm castors





AX 190
Imprest Trolley
 1270 x 630 x 1330mm
 Powder coated finish
 4 Shelves - mesh sides and back
 Full bumper surround / 200mm castors
 Available with 6 shelves



AX 191
Imprest Trolley
 1270 x 630 x 1330mm
 Powder coated finish
 4 Shelves - enclosed sides and back
 Full bumper surround / 200mm castors
 Available with 6 shelves

Sizes available
 R2500 - 16 compartment
 R2510 - 18 compartment
 R2520 - 24 compartment
 R2530 - 32 compartment
 (pictured)



R2530
Residents Clothing Distribution Trolley
 Multi compartment trolley designed for distribution of personal linen for residents. Features double sided access. Comprises shelf ledges, vertical rods, hanging rods and push handle.



AX 723
Personal Distribution Trolley
 960 x 680 x 1310mm
 Powder coated frame
 20 bays / 20 individual tubs included
 Additional sizes available



AX 199
Removable Basket Trolley
 750 x 490 x 900mm
 Nylon coated baskets - 200mm deep
 Push handle / 125mm castors
 BASK 0202 Replacement Basket



AX 205
Fixed Basket Trolley
 750 x 490 x 900mm
 Nylon coated baskets - 200mm deep
 Push handle / 125mm castors



AX 147
Clearing Trolley - 2 dished shelves
 900 x 490 x 900mm
 Stainless steel construction



AX 625
Clearing Trolley
 900 x 490 x 900mm
 Stainless steel construction
 2 break shelves



AX 662
Light Weight Tray Trolley
 680 x 550 x 1500mm - 2 x 10 capacity
 Tubular stainless steel construction
 AX 658 - 1 x 10 capacity
 AX 665 - 3 x 10 capacity



AX 150
Clearing Trolley - 3 dished shelves
 900 x 490 x 900mm
 Stainless steel construction



AX 627
Clearing Trolley
 900 x 490 x 900mm
 Stainless steel construction
 3 break shelves



AX 699
Light Weight Tray Trolley
 815 x 480 x 1500mm - 2 x 10 capacity
 Depending on tray size
 Stainless steel construction



Dished shelf
 5mm deep




Break shelf
 2mm deep



AX 654
Mobile Preparation Table
 1200 x 760 x 900mm
 Stainless construction
 125mm braking castors

Other size available
 AX 651 - 1500 x 760 x 900mm
 AX 648 - 1800 x 760 x 900mm
 Static Units
 AX 645 - 1200 x 760 x 900mm
 AX 642 - 1500 x 760 x 900mm
 AX 639 - 1800 x 760 x 900mm

CONSULTING ROOMS



Our consulting room range offers an extensive selection of product options to accommodate all types of consulting rooms as well as coordinating with the surrounding décor. The inclusion of fully electrically operated examination couches addresses patient safety for ingress and egress to the couch and allows flexibility for the physician to treat the patient at a comfortable height.

The options of stainless steel or colour coordinated powder coated equipment is available and accommodates different price levels.



AX3610

OmegaE Electric 3 Section Couch

- All electric operation

- Electric operated variable height
- Electric operated head and foot section
- Height range: 490 - 900mm
- Head section: positive and negative tilt
- Size: 1900mm long x 700mm wide
- Padded, upholstered top, with face cut out
- 125mm castors - fully lockable
- Safe working load 200kg



AX 3600

OmegaE Electric 2 Section Couch

- All electric operation

- Electric operated variable height
- Electric operated head and foot section
- Height range: 490 - 900mm
- Head section: positive and negative tilt
- Size: 1900mm long x 700mm wide
- Padded, upholstered top, with face cut out
- 125mm castors - fully lockable
- Safe working load 200kg
- Optional electric foot control



SEERS SX 2007

Sterling 2 Section Couch

- Electric variable height with hand control as standard
- Gas strut assisted adjustable backrest from -25° to +75° positive
- Height range from 450 - 970mm
- Standard width 700mm
- Strong and rigid frame design
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- 100mm castors system with locking brakes
- Safe working load 200kg



SEERS SX 3008

Sterling 3 Section Couch

- Electric variable height with hand control as standard
- Gas strut assisted adjustable headrest from -25° to +40° positive
- Gas strut assisted adjustable footrest from 0° to +70° positive
- Height range from 450 - 970mm
- Standard width 700mm
- Strong and rigid frame design
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- 100mm castors system with locking brakes
- Safe working load 200kg

CONSULTING ROOMS



AX 357G
Examination Couch
 1810 x 600 x 820mm
 Adjustable head section
 Tubular frame
 EXAM 0354 - Square Tube



AX 363
Sliding Arm Board
 To suit suit all models of
 Examination Couches



AX 381
Paper Roll Holder
 To suit suit all models of
 Examination Couches



AX 385
Examination Couch Heavy Duty
 1810 x 760 x 820mm
 Adjustable head section
 Reinforced frame structure
 SWL 250kg

AX 382
3 Drawer Unit (optional)
 Modular drawer unit to
 suit couch AX 385
 Drawers on metal roller
 runners

Examination Couch and Drawer Unit combinations

- AX 385 DLHS
Left Side Drawer Unit
- AX 385 DRHS
Right Side Drawer Unit

PICTURED WITH OPTIONAL DRAWER UNIT



AX 375G
One Step Stool
 Step height 230mm
 Epoxy coated - non slip tread



AX 393G
Two Step Stool
 Step height 230mm
 Epoxy coated - non slip tread



AX 378
One Step Stool with handle rail
 Step height 250mm
 Epoxy coated - non slip tread
 Handle to one end



AX 402RHR - Right
Two Step Stool with hand rail
 Step heights 200 / 400mm
 Epoxy coated - non slip tread
 Hand rail available left or right
 AX 402LHR - Left

POWDER COATED



AX 321
Instrument Trolley
490 x 490 x 900mm



AX 324
Instrument Trolley - 1 Drawer
490 x 490 x 900mm



AX 327
Instrument Trolley - 2 Drawer
490 x 490 x 900mm



AX 333
Instrument Trolley - 1 Drawer
750 x 490 x 900mm



AX 336
Instrument Trolley - 2 Drawer
750 x 490 x 900mm

NOVA



AX 342
Instrument Trolley
490 x 490 x 900mm
Powder coated finish -
Stainless steel shelves



AX 345
Instrument Trolley - 1 Drawer
490 x 490 x 900mm
Powder coated finish
Stainless steel shelves



AX 348
Instrument Trolley - 2 Drawer
490 x 490 x 900mm
Powder coated finish
Stainless steel shelves



AX 349
Instrument Trolley - 3 Drawer
490 x 490 x 900mm
Powder coated finish
Stainless steel shelves

STAINLESS STEEL



AX 033
Instrument Trolley
490 x 490 x 900mm
Stainless steel construction
75mm castors



AX 054
Instrument Trolley - 1 Drawer
490 x 490 x 900mm
Stainless steel construction
75mm castors



AX 066
Instrument Trolley - 2 Drawer
490 x 490 x 900mm
Stainless steel construction
75mm castors



AX 076
Instrument Trolley - 3 Drawer
490 x 490 x 900mm
Stainless steel construction
125mm castors

CONSULTING ROOMS



AX 468
Blood Collection Chair
800 x 1050 x 750mm
Powder coated finish



AX 473
Blood Collection Chair
Stainless steel Sliding Footrest
800 x 1050 x 750mm
Powder coated finish



AX 411
Privacy Screen 3 fold
Panel size 900 x 1400mm
Unfolded 2700 x 1700mm
Powder coated frame
Includes curtains
AX 414CURT
Replacement curtain / 3 fold
AX 423CURT
Replacement curtain / 2 fold



AX 262
Surgeon Stool Hand Operated
Height 530 - 710mm
No backrest



AX 265
Surgeon Stool Hand Operated
Height 530 - 710mm
Backrest



AX 2710
Saddle Stool
Hand operated height adjustment
Height range 530 - 730mm



AX 420
Privacy Screen 2 fold
Panel size 900 x 1400mm
Unfolded 1800 x 1700mm
Powder coated frame
Includes curtains



AX 089
Tub Trolley
490 x 490 x 900mm
4 x Removable tubs
Laminate top
Lift-up side shelf



AX 2382
Examination Light Stand
Multi Purpose
Powder coated finish
5 Star base - 50mm castors



AX 9122
Tilt Bins - 6 Units



AX 9124
Tilt Bins - 5 Units



AX 9126
Tilt Bins - 4 Units



AX 9128
Tilt Bins - 3 Units



AX 9144
Document Holder
260 x 254 x 102mm



AX 9142
Document Holder
254 x 267 x 51mm




AX 9140
Document Holder
445 x 165 x 76mm



AX 9136
Document Holder
368 x 165 x 76mm
Set of 3



AX 9146
Document Holder
375 x 178 x 102mm



Provita Series 1
Clean, modern design, Provita examination lamps have the latest LED technology for economic long term operation.

LED bulbs and the specifically developed reflectors focus clean light where you need it. Optimal heat dissipation utilising a double wall lamp head housing reduces any heat output.

Easy and safe positioning of lamp head with integrated ergonomic handle.




PROVL 100011A
LED Examination Lamp with Gooseneck Arm
Working radius: approx 800mm
Bulb life expectancy: 20,000 hours
Switch at lamp head
Total weight: 2.6kg



PROVL 100031A
LED Examination Lamp with Joint Articulated Arm
Working radius: approx 800mm
Bulb life expectancy: 20,000 hours
Switch at lamp head
Total weight: 2.8kg



PROVL 100081A
LED Examination Lamp with Double-Joint Articulated Arm
Working radius: approx 1100mm
Bulb life expectancy: 20,000 hours
Switch at lamp head
Total weight: 2.8kg



AX 2382
Mobile Examination Lamp with Fixed Height Stand
(lamp not included)

Stable tilt resistant roller base
Base diameter 635mm
Castors 50mm with brake
Steel construction
White powder coat finish

PRODUCT GUIDE

PRODUCT	CODE	PAGE
Airway Trolley	AX 106	13
Anaesthetic Trolley - four drawer	AX 102	05
Anaesthetic Trolley - four drawer	AX 104	05
Anaesthetic Trolley - six drawer	AX 105	05
Arm Table - adjustable height	AX 113	13
Arm Table - fixed height	AX 111	13
AV Equipment Trolley - modular	AX 841	15
Babies Bath - adjustable height	AX 610	25
Basket Trolley - fixed	AX 205	30
Basket Trolley - removable	AX 199	30
Bassinet Trolley - single adjustable height	AX 609	25
Bassinet Trolley - single	AX 582	25
Bassinet Trolley - twin	AX 597	25
Bassinet Trolley - single with cabinet	AX 612	25
Bed Cradle	AX 501	19
Blood Collection Chair	AX 468	36
Blood Collection Chair - sliding footrest	AX 473	36
Bowl Stand - single	AX 211	15
Bowl Stand - single SML	AX 248	15
Bowl Stand - double	AX 214	15
Bulk Collection Trolley	AX 634	29
Charts Trolley - 12 charts	AX 990	19
Charts Trolley - 24 charts	AX 993	19
Charts Trolley - 36 charts	AX 994	19
Chart Holder A4	AX 996	19
Clearing Trolley - 2 break shelves	AX 625	31
Clearing Trolley - 2 dished shelves	AX 147	31
Clearing Trolley - 3 break shelves	AX 627	31
Clearing Trolley - 3 dished shelves	AX 150	31
CSSD Case Cart	AX 692	11
CSSD Kimguard Drying Rack	AX 0681	11
CSSD List Cart	AX 278	11
CSSD Packing Table	AX 0680	11
CSSD Preparation Table - adjustable height	AX 0688	11
CSSD Transportation Trolley	AX 0684	11
CSSD Work Station	AX 0666	10
CSSD Work Station Accessories		
Accessory Shelf	AX 0672	10
CPU Holder	AX 0686	10
Keyboard Shelf	AX 0687	10
Lower Storage Shelf	AX 0689	10
Power Board Outlet 6 plug	AX 0690	10
Sterile Wrap Rack	AX 0685	10
Steri-Peel Holder	AX 0691	10
Distribution Trolley - Clothing	R2530	30
Distribution Trolley - Personal	AX 723	30
Document Holder - 260 x 254 x 102mm	AX 9144	23
Document Holder - 254 x 267 x 51mm	AX 9142	23
Document Holder - 445 x 165 x 76mm	AX 9140	23
Document Holder / set of 3 - 368 x 165 x 76mm	AX 9136	23
Document Holder - 375 x 178 x 102mm	AX 9146	23
Files Trolley	AX 655	19
Emergency Treatment Trolley	AX 115	06
Emergency Treatment Trolley	AX 117	06
Emergency Treatment Trolley	AX 119	06
Equipment Trolley - Modular	AX 840	15
Examination Couch	AX 357G	34
Examination Couch Heavy Duty	AX 385	34
Examination Couch Accessories		
3 Drawer Unit	AX 382	34
Paper Roll Holder	AX 381	34
Sliding Arm Board	AX 363	34
Examination Couch - Electric 2 section	AX 3600	33
Examination Couch - Electric 3 section	AX 3610	32
Examination Couch - Electric 2 section	SEERS SX 2007	33
Examination Couch - Electric 3 section	SEERS SX 3008	33
Garbag Trolley	AX 519	18
Garbag Wall Mount Bracket	AX 522	19
Garment Rack	AX 175	18
Gas Bottle Trolley - C size single	AX 537	18
Gas Bottle Trolley - C size single stainless	AX 540	18
Gas Bottle Trolley - D size single	AX 546	18
Gas Bottle Trolley - D size single stainless	AX 549	18
Gas Bottle Trolley - G size single	AX 528	18

PRODUCT	CODE	PAGE
Gas Bottle Trolley - G size single stainless	AX 531	18
Gas Bottle Holder - hook over	AX 534	19
Gas Bottle Holder - hook over	AX 543	19
Gas Bottle Trolley - stainless steel 9 x C size	AX 524	18
Gas Bottle Trolley - stainless steel 12 x C size	AX 525	18
Gas Bottle Trolley - stainless steel multi D/C	AX 524SS	18
Gas Bottle Trolley - stainless steel C size	AX 525SS	18
Imprest Trolley	AX 191	30
Imprest Trolley - mesh	AX 190	30
Instrument Trolley - no rail	AX 001	02
Instrument Trolley - no rail	AX 003	02
Instrument Trolley - no rail	AX 006	02
Instrument Trolley - no rail	AX 009	02
Instrument Trolley - rail 3 sides	AX 033	02
Instrument Trolley - rail 3 sides	AX 036	02
Instrument Trolley - rail 3 sides	AX 039	02
Instrument Trolley - rail 3 sides	AX 042	02
Instrument Trolley - one drawer	AX 054	02
Instrument Trolley - one drawer	AX 057	02
Instrument Trolley - one drawer	AX 060	02
Instrument Trolley - one drawer	AX 063	02
Instrument Trolley - two drawer	AX 066	02
Instrument Trolley - two drawer	AX 069	02
Instrument Trolley - two drawer	AX 072	02
Instrument Trolley - two drawer	AX 075	02
Instrument Trolley - two half drawers	AX 073	03
Instrument Trolley - three drawer	AX 076	03
Instrument Trolley - four drawer	AX 080	03
Instrument Trolley - multi drawer	AX 078	05
Instrument Trolley - Powder coated / stainless drawers		
Instrument Trolley	AX 342	35
Instrument Trolley - one drawer	AX 345	35
Instrument Trolley - two drawer	AX 348	35
Instrument Trolley - three drawer	AX 349	35
Instrument Trolley - Powder Coated		
Instrument Trolley	AX 321	35
Instrument Trolley - one drawer	AX 324	35
Instrument Trolley - two drawer	AX 327	35
Instrument Trolley - one drawer	AX 333	35
Instrument Trolley - two drawer	AX 336	35
ICU Work Station SML	AX 701	15
ICU Work Station LRG	AX 702	15
IV Set Up Trolley	AX 145	06
IV Stand	AX 217	17
IV Stand - heavy base	AX 219	17
IV Stand - polymer base	AX 216	17
IV Stand - POGO heavy base	AX 215	17
IV Stand - transfusion pump	AX 220	17
IV Stand - smart IV	AX 2380	17
IV Stand Accessories		
Handle and Clamp	AX 233	17
Storage Basket and Clamp	AX 231	17
Vertical Oxy Holder and Clamp	AX 232	17
Kick Bucket Stand - combined	AX 240	15
Kick Bucket Stand	AX 241	15
Kick Bucket - 11.4 litre	AX 244	15
Kick Bucket Lid	AX 247	15
Lamp Stand - mobile	AX 2382	37
Lamp Stand - multi purpose	AX 2382	36
Laundry Trolley - Moist 320 litre	AX 729	29
- Spring loaded base only	AX 731	29
Laundry Trolley - moist 450 litre	AX 733	29
Laundry Trolley - wet and dry	AX 726	29
Lead Apron Rack	AX 177	18
LED Examination Lamp - Double-Joint Arm	PROVL100081A	37
LED Examination Lamp - Gooseneck Arm	PROVL100011A	37
LED Examination Lamp - Joint Articulated Arm	PROVL100031A	37
Linen Skip - Double	AX 294	27
Linen Skip - Double Foot Operated Lid	AX 300	27
Linen Skip - Half Height Shelf	AX 320	27
Linen Skip - Single	AX 286	27
Linen Skip - Single Foot Operated Lid	AX 291	27
Linen Skip - Triple	AX 303	27
Linen Skip - Triple Foot Operated Lid	AX 307	27

PRODUCT	CODE	PAGE
Linen Skip Powder Coated		
Linen Skip - single foot operated lid	AX 0291 PC	28
Linen Skip - double foot operated lid	AX 0300 PC	28
Linen Skip - triple foot operated lid	AX 0307 PC	28
Linen Skip Options		
Soft Close Lid System	AX 319	27
Swivel Hooks	AX 315	27
Linen Trolley - soiled	AX 0730	29
Lithotomy Trolley	AX 260	13
Mayo Table - slide base	AX 253	13
Mayo Table - two leg base	AX 252	13
Mayo Table - three leg base	AX 254	13
Mayo Table - four leg base	AX 256	13
Medication Trolley - one drawer	AX 677	07
Medication Trolley - three drawer	AX 676	07
Mortuary / Concealment Trolley	AX 704	22
Mortuary / Concealment Trolley Cover	AX 707	22
Nesting Tables - set of three	AX 021	13
Nesting Table - large	AX 024	13
Nesting Table - medium	AX 027	13
Nesting Table - small	AX 030	13
Pan & Bottle Rack - 4	AX 554	19
Pan & Bottle Rack - 6	AX 552	19
Patient Trolley - fixed height	AX 708	22
Plaster Trolley - two drawer	AX 4990	13
Preparation Table - mobile	AX 654	31
Privacy Screen 2 fold	AX 420	36
Privacy Screen 3 fold	AX 411	36
Privacy Screen replacement curtain / 2 fold	AX 423CURT	36
Privacy Screen replacement curtain / 3 fold	AX 414CURT	36
Shower Chair	AX 432	20
Shower Chair - swing back arms	AX 432SA	20
Shower Chair - fixed rear wheels	AX 434	21
Shower Chair - heavy duty	AX 433	21
Shower Chair - bariatric	AX 435	21
Shower Chair Accessories		
IV Pole and Clamp	AX 460	21
Leg Extension	AX 447	21
Safety Arms - Wrap Around	AX 455	21
Vertical Oxy Holder and Clamp	AX 462	21
Spectra 2 Shelf Trolley	AX S0001	08
Spectra 1 Drawer Trolley	AX S0010	08
Spectra 2 Drawer Trolley	AX S0020	08
Spectra Multi Drawer Trolley	AX S0030	08
Step Stool-one step epoxy coated	AX 375G	34
Step Stool-two step epoxy coated	AX 393G	34
Step Stool-one step epoxy coated/handle rail	AX 378	34
Step Stool-two step epoxy coated/handle rail right	AX 402RHR	34
Step Stool-two step epoxy coated/handle rail left	AX 402LHR	34
Step Stool-one step stainless/aluminium tread	AX 387	14
Step Stool-one step stainless/aluminium tread	AX 390	14
Step Stool-two step stainless/aluminium tread	AX 405	14
Step Stool-one step stainless/alum tread/hand rail	AX 391	14
Storage Trolley	AX 094	29
Storage Unit Mobile - 3 trays	AX M2500	09
Storage Unit Mobile - 5 trays	AX M1400	09
Storage Unit Mobile - 7 trays	AX M8000	09
Storage Unit Mobile - 10 trays	AX M2000	09
Storage Unit Static - 7 trays	AX M8500	09
Surgeons Stool foot operated	AX 267	14
Surgeons Stool foot operated	AX 270	14
Surgeons Stool hand operated	AX 262	14
Surgeons Stool hand operated	AX 265	14
Surgeons Stool - hydraulic	AX 273	14
Surgeon Stool - saddle type	AX 2710	14
Resuscitation Trolley - 4 drawers	AX 107	03
Resuscitation Trolley - standard 5 drawers	AX 108	05
Resuscitation Trolley - 6 drawers	AX 110	05
Theatre Stool	AX 277	14
Theatre Supply Storage Cart	AX 636	15
Tilt Bin - 3 Bin Unit	AX 9128	23
Tilt Bin - 4 Bin Unit	AX 9126	23
Tilt Bin - 5 Bin Unit	AX 9124	23
Tilt Bin - 6 Bin Unit	AX 9122	23

PRODUCT	CODE	PAGE
Tray Trolley	AX 699	31
Tray Trolley - Tubular	AX 662	31
Trays - Sixty Forty System		
50mm deep tray	AX M0601	08
Interlocking divider 600mm	AX M1661	08
Interlocking divider 400mm	AX M1641	08
100mm deep tray	AX M0605	08
Interlocking divider 600mm	AX M1665	08
Interlocking divider 400mm	AX M1645	08
200mm deep tray	AX M0602	08
Interlocking divider 600mm	AX M1662	08
Interlocking divider 400mm	AX M1642	08
Stainless Steel Trolley Accessories		
Drawer insert containers	AX 120	07
Drop down work surface	AX 126	07
IV pole and clamp	AX 129	07
Power board	AX 135	07
Push handle	AX 143	07
Storage panel	AX 132	07
Self closing drawer runner	AX 738	07
Vertical oxy bottle holder and clamp	AX 123	07
Waste bin and bracket	AX 142	07
Tub Trolley	AX 089	36
Urinal Bottle Holder - nylon coated	AX 510	19
Urinal Bottle Holder - stainless steel	AX 513	19
Utility Trolley	AX 091	29

GM Medical Pty Ltd
28-30 Hightech Place
Lilydale VIC 3140
T 03 9735 0111
F 03 9735 2104
gmmedical.com.au

ParagonCare

Paragon Care Group of Companies

